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DNA isolation

Ethical approval was granted by...(if needed)

Blood / Tissue: Genomic DNA was extracted from ... using DNeasy Blood & Tissue Kit (QIAGEN)

FFPE: Genomic DNA was extracted from ... using QIAamp DNA FFPE tissue kit (QIAGEN)

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Library preparations and Sequencing

The genomic DNA was sent to **Shanghai Biotechnology Corporation** for the de novo wholegenome sequencing. Library was prepared according to the TruSeq DNA Sample Preparation Guide (Illumina, 15026486 Rev.C) . The quantity of libraries were assessed by Qubit® 2.0 Fluoromete. The quality and size of libraries were measured by 2100 Bioanalyzer High Sensitivity

DNA Assay according to the reagent kit guide. For Illumina sequencing, the qualified libraries were applied to 2×150 bp paired-end sequencing on Illumina HiSeq X-ten platform (Illumina).

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Reads Mapping

FASTQ files were aligned to human reference genome (hg19/GRCh37) by BWA v0.7.13. The aligned files (sam/bam format files) were sorted by samtools (1.3) firstly, then duplicates were flagged by using Picard (v2.2.4). By using GATK v3.5, reads were locally realigned and base qualities were recalibrated. Finally, mapping statistics include coverage and depth were generated from recalibrated files by BEDTools (v2.16.1) and in-house perl/python scripts.

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Variant calling and annotation

Variants (SNVs and indels) were genotyped from recalibrated BAM files using the multi-sample processing mode of the Unified Genotyper tool from the GATK. Then VQSR (Variant Quality Score Recalibration) was used to reduce false positive of variant calling. Copy Number Variation and Structural Variation were identified by CNVnator v0.2.7 and LUMPY v0.2.13, respectively. SNVs and indels were annotated using ANNOVAR software against multiple databases including HGVS variant description, population frequency, disease or phenotype and variant functional prediction.